

Sterilization & Drying

Electric Ovens Dry Heat Sterilizers

Panasonic sterilizers provide constant high-temperature conditions for dry heat sterilization. Panasonic electric ovens combine an easy-to-use design and strict safety features to facilitate a host of applications.




Panasonic MOV Series

Electric Ovens
and Dry Heat Sterilizers

Accurate, High-Temperature Equipment for the New Generation of Scientific Research.

Panasonic has always aimed to provide research support equipment that offers complete satisfaction to its users. Inspired by the search for even higher precision and greater flexibility of control, Panasonic presents the new MOV Series.



Search to the future

Microprocessor PID temperature control system guarantees accurate temperature environment

The microprocessor PID (Proportional, Integrated and Differential) temperature control system ensures accurate inside temperature.

With less offset or overshoot, exact control is possible. This system provides the high-temperature environment that exactly meets experimentation requirements.



Forced air circulation system ensures stable temperatures accurate to within ±4.0°C


Fan circulation ensures that deviations in cabinet temperature are kept within ±4.0°C (at 200°C).

The MOV Series can be widely used for basic to applied experimentation in the areas of scientific, industrial and environmental testing.

Sheathed heater ensures durability and safety

A sheathed heater is incorporated in the heater section. The heating element is wrapped in a magnesium oxide insulating material and covered with metal protection tube. With conventional wire heating elements (bare heaters), gases or dust can cause corrosion, resulting in loss of heating capacity and electrical leakage.

With its durability and high chemical resistance, the sheathed heater ensures safer, more stable operation without the risk of electrical leakage.



Attractive design

Panasonic believes that laboratory equipment should be attractive as well as functional. The MOV Series features a future-oriented design, with rounded corners, door handles that blend with the main body, and a flat control panel.



User-oriented design for easy operation

The control panel has soft-touch keys and bright, green digital LED display that allows easy confirmation of temperature and remaining operation time. Other advantages of the design include a soft-latched door handle integrated with the door, an observation window for checking conditions inside the cabinet, two exhaust vents (shared with an access vent) on the top of the unit, and a stainless-steel (SUS-304) interior to guarantee durability and superior resistance to chemicals.

Alarm and safety functions

A comprehensive range of alarm and safety devices is included as standard in the MOV Series. And the built-in alarm jack allows remote control of the unit.



Alarm and safety functions	Trigger	Alarm notification method and operation
Automatic set temperature alarm	More than 10°C deviation from set temperature.	Lamp, LED, buzzer.
Independent overheating protection circuit	Abnormal temperature increase above upper limit.	While not in operation: buzzer. While in operation: lamp, LED and buzzer. Independent circuit switches off heater and fan motor.
Overtemp. safety system for control section	Ambient temperature of main part (base) of control section exceeds 65°C.	Lamp, LED, buzzer. Heater and fan motor switched off.
Keylock switch	Keylock switch on.	Key input not possible (excl. buzzer and call key).
Memory backup function	During power failure, breaker cut or when unit is not connected to power supply.	Stores operation patterns for set temperatures/times.

Malfunction Monitor (Self diagnosis function)

Should a malfunction occur, it is diagnosed and indications are given on the digital display.

Indication	Meaning
E 01	Temperature sensor open
E 02	Triac open
E 03	Triac shortcircuit
E 04	Relay shortcircuit
E 05	Relay open, fan motor shortcircuit, heater shortcircuit, independent over-heating protection circuit activated and overtemp. safety system is ON.

Electric Ovens

MOV-112F/212F
MOV-112/212

Future-oriented Ovens That Look Good and Improve Efficiency.

Four models that feature natural convection and forced air circulation systems to create environments for a wide variety of experiments. Designed for ease of use and safety.

Microprocessor timer function

Panasonic has included a microprocessor timer function, so operation times can be set up to a maximum of 99 hours and 59 minutes. The combination of auto start and auto stop provides operating patterns suited to a wide variety of applications. The auto stop operates the timer when the heater is on, or when the set temperature has been reached. A buzzer indicates the end of timer operation.

Forced air circulation system (MOV-112F/212F)

Sirocco fan circulation keeps variations in inner cabinet temperature within $\pm 4^{\circ}\text{C}$ at 200°C . Compared with natural convection, quicker drying is possible. And Panasonic's unique fan motor not only circulates hot air in the cabinet but also keeps the motor cool, improving the reliability and safety of the motor.

Natural convection system (MOV-112/212)

Natural convection is best for drying very small samples and fine particles which would be scattered by a fan. This system can be used for high-temperature applications up to 250°C .

Fast heat treatment possible at an even temperature.

Ambient temperature $+5^{\circ}\text{C}$ to 200°C
(Forced air circulation system)

Effective capacity

90 liters

MOV-112F



Effective capacity

150 liters

MOV-212F

Versatile operating temperature range without damaging samples.

Ambient temperature $+5^{\circ}\text{C}$ to 250°C
(Natural convection system)

Effective capacity

97 liters

MOV-112




Effective capacity

157 liters

MOV-212

Dimensional Figures

MOV-112F




Unit: mm

Performance Data

MOV-112F


Pull-up characteristics

No load, exhaust vent closed.



Dimensional Figures

MOV-112




Unit: mm

Performance Data


MOV-112

Pull-up characteristics

No load, exhaust vent closed.



MOV-212F




Unit: mm


MOV-212F

Pull-up characteristics

No load, exhaust vent closed.



MOV-212




Unit: mm

MOV-212

Pull-up characteristics

No load, exhaust vent closed.



Dry Heat Sterilizers

MOV-112S/212S

Constant Temperature Environments for Dry Heat Sterilizing and Efficient Laboratory Work.

MOV Series models provide many advantages:
PID precision temperature control is adjustable to within $\pm 1^{\circ}\text{C}$; Forced air circulation keeps cabinet temperatures even to within $\pm 4^{\circ}\text{C}$; The new microprocessor timer helps correct sterilizing time.

The microprocessor PID control system guarantees accurate sterilizing temperatures.

The built-in sheathed heater offers superior durability and safety.

Forced air circulation with quick pull-up and precise temperature distribution

The time needed to reach the set temperature is shorter with a forced air circulation system compared with natural convection. The temperature is even throughout the cabinet—deviations are kept within $\pm 4^{\circ}\text{C}$ at 200°C . And the fan motor is specially designed to circulate the air in the cabinet and prevent the motor from overheating at the same time.




Built-in microprocessor timer guarantees accurate sterilization time

The built-in microprocessor timer (max. setting 99 hours and 59 minutes) accurately adjusts sterilization time. A consecutive display of the inner cabinet temperature is provided for quick confirmation of the time when samples should be removed. And a buzzer sounds at the end of timer operation. In combination with auto start and auto stop, operation patterns can be set for a wide range of applications.

Standard dry heat sterilizing times (as indicated by the Japanese Pharmacopoeia)

$135^{\circ}\text{C} - 145^{\circ}\text{C}$	3 – 5 hours
$160^{\circ}\text{C} - 170^{\circ}\text{C}$	2 – 4 hours
$180^{\circ}\text{C} - 200^{\circ}\text{C}$	0.5 – 1 hours

Microprocessor timer setting patterns:



Accurate temperature and time — essential for dry heat sterilization.
Ambient temperature $+5^{\circ}\text{C}$ to 200°C




Effective capacity
90 liters
MOV-112S

Effective capacity
150 liters
MOV-212S


Dimensional Figures

MOV-112S



Unit: mm

MOV-212S




Unit: mm

Performance Data


MOV-112S Pull-up characteristics

No load, exhaust vent closed.



MOV-212S Pull-up characteristics

No load, exhaust vent closed.



Specifications

Type	Electric Ovens				Dry Heat Sterilizers		
Model No.	110V, 60Hz	MOV-112F-PT	MOV-212F-PT	MOV-112-PT	MOV-212-PT	—	MOV-212S-PT
	220V, 50/60Hz	MOV-112F-PK	MOV-212F-PK	MOV-112-PK	MOV-212-PK	MOV-112S-PK	MOV-212S-PK
	230V, 50Hz (CE)	MOV-112F-PE	MOV-212F-PE	MOV-112-PE	MOV-212-PE	MOV-112S-PE	MOV-212S-PE
	240V, 50Hz	MOV-112F-PU	MOV-212F-PU	MOV-112-PU	MOV-212-PU	MOV-112S-PU	MOV-212S-PU
External dimensions (W x D x H)*1	580 x 595 x 820 (mm) 22.8 x 23.4 x 32.3 (inch)	730 x 645 x 870 (mm) 28.7 x 25.4 x 34.3 (inch)	580 x 595 x 820 (mm) 22.8 x 23.4 x 32.3 (inch)	730 x 645 x 870 (mm) 28.7 x 25.4 x 34.3 (inch)	580 x 595 x 820 (mm) 22.8 x 23.4 x 32.3 (inch)	730 x 645 x 870 (mm) 28.7 x 25.4 x 34.3 (inch)	
Internal dimensions (W x D x H)	450 x 450 x 450 (mm) 17.7 x 17.7 x 17.7 (inch)	600 x 500 x 500 (mm) 23.6 x 19.7 x 19.7 (inch)	450 x 475 x 450 (mm) 17.7 x 18.7 x 17.7 (inch)	600 x 525 x 500 (mm) 23.6 x 20.7 x 19.7 (inch)	450 x 450 x 450 (mm) 17.7 x 17.7 x 17.7 (inch)	600 x 500 x 500 (mm) 23.6 x 19.7 x 19.7 (inch)	
Effective capacity	90 liters (3.2 cu.ft.)	150 liters (5.3 cu.ft.)	97 liters (3.4 cu.ft.)	157 liters (5.5 cu.ft.)	90 liters (3.2 cu.ft.)	150 liters (5.3 cu.ft.)	
Exterior finish	Painted steel						
Interior finish	Stainless steel plate						
Insulation	Glass wool						
See-through window	Reinforced triple-glass window (t = 5 mm)						
Shelves	Stainless-steel plate, stainless-steel wire (adjustable)						
	2	3	2	3	2	3	
Air exhaust vent	Two on top plate (23 mm inside dia.)						
Heating system	Forced air circulation system		Natural convection system		Forced air circulation system		
Temperature control system	Microprocessor PID control						
Sensor	Thermo couple						
Temperature setting	Digital setting (adjustable range: ±1°C)						
Timer	Auto start, Auto stop 00:00 – 99:59						
Temperature/Timer display	Digital LED display						
Heater (Sheathed heater)	1.1kW	1.2kW	1.1kW	1.3kW	1.1kW	1.2kW	
Interior fan	Sirocco fan dia. 149 mm		—		Sirocco fan dia. 149 mm		
Exterior fan	Propeller fan 107 mm		—		Propeller fan 107 mm		
Power source	50/60Hz, cord approx. 2m						
Max. power consumption	Approx. 1.1kW	Approx. 1.2kW	Approx. 1.1kW	Approx. 1.3kW	Approx. 1.1kW	Approx. 1.2kW	
Temperature range	Ambient temperature +5°C to 200°C		Ambient temperature +5°C to 250°C		Ambient temperature +5°C to 200°C		
Temperature controllability	±1 deg.				—		
Temperature uniformity	±4°C (at 200°C)		±10°C (at 200°C)		±4°C (at 200°C)		
Weight	50kg (110.2 lbs.)	66kg (145.5 lbs.)	47kg (103.6 lbs.)	63kg (138.9 lbs.)	50kg (110.2 lbs.)	66kg (145.5 lbs.)	
Alarm and safety function	Overcurrent breaker, automatic set temperature alarm (set point +10°C), independent overheat protection circuit, overtemperature safety system at control section (triggered at 65°C), self diagnosis.						

Appearance and specifications are subject to change without notice.

When the setting value is below 100°C, large deviations in cabinet temperature will occur.

*1 Exterior dimensions of main cabinet only - see dimension drawings showing handles and other external projections.



Panasonic Healthcare Co., Ltd., Biomedical Business Unit is certified for:
Quality management system: ISO9001
Medical devices quality management system: ISO13485



Panasonic Healthcare Co., Ltd., Biomedical Business Unit is certified for:
Environmental management system: ISO14001

DISTRIBUTED BY:

Panasonic®

<http://panasonic.net/healthcare/biomedical/>

Printed in Japan 041-2013-03